
Der **Kombinations-Test** neuer Coronaviren (SARS-CoV-2-Antigene), Influenza A-Virus-Antigene, Influenza B-Virus-Antigene und RSV-Antigene

4 in 1

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WHAT DOES THE KIT TEST?

The fluorecare® Combination Test **(4 in 1)** for SARS-CoV-2 and Influenza A/B virus and RSV antigens can be used for the simultaneous qualitative detection and differentiation of novel coronaviruses (SARS-CoV-2 antigens), influenza A virus antigens, influenza B virus antigens and RSV antigens in nasal swab specimens in vitro.

This kit finds its use in the diagnoses of infection caused by SARS-CoV-2 coronavirus (COVID-19), as well as infection caused by influenza A or B or RSV virus.

For in vitro diagnostic use only. For 1-pack and 5-pack only for self-testing, 25-pack only for professional use.

Der **Kombinations-Test** neuer Coronaviren (SARS-CoV-2-Antigene), Influenza A-Virus-Antigene, Influenza B-Virus-Antigene und RSV-Antigene

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ABOUT COVID-19

The new coronaviruses belong to the β -genus. COVID-19 is an acute respiratory infectious disease. All people are susceptible to it. Currently, people infected with the new coronavirus are the main source of infection. Asymptomatic infected individuals may also be a source of infection. Based on current epidemiologic studies, the incubation period of new coronavirus is 1 to 14 days, most commonly 3 to 7 days. The most common manifestations are fever, fatigue, and irritable cough. In a few cases, nasal congestion, runny nose, sore throat, muscle aches, and diarrhea occur.

Der **Kombinations-Test** neuer Coronaviren (SARS-CoV-2-Antigene), Influenza A-Virus-Antigene, Influenza B-Virus-Antigene und RSV-Antigene

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TEST PRINCIPLE

The SARS-CoV-2, influenza A, influenza B, and RSV antigens are qualitatively detected in nasal swab specimens of the subjects by the colloidal gold method. After the sample is added, the SARS-CoV-2 antigen (or influenza A/B and RSV) in the sample to be tested is combined with the SARS-CoV-2 antigen (or influenza A/B and RSV) antibody labeled with colloidal gold on the binding pad to form the SARS-CoV-2 antigen (or influenza A/B and RSV) antibody-colloidal gold complex. Because of chromatography, the SARS-CoV-2 antigen (or influenza A/B and RSV) antibody-colloidal gold complex diffuses along the nitrocellulose membrane. Within the detection line region, the SARS-CoV-2 antigen (or influenza A/B and RSV) -antibody complex binds to the antibody in the detection line region and shows a purple-red band. Colloidal gold labeled SARS-CoV-2 antigen (or influenza A/B and RSV) diffuses to the quality control line region (C) and is captured by sheep anti-mouse IgG and shows red bands. When the reaction is complete, you can interpret the results by visual observation.

Der **Kombinations-Test** neuer Coronaviren (SARS-CoV-2-Antigene), Influenza A-Virus-Antigene, Influenza B-Virus-Antigene und RSV-Antigene

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1s test

Product ID # **MF-71-1**

- 1 Test cassette
- 1 Sterile nasal swab
- 1 Prefilled sample Treatment Tube
- 1 Instruction manual D/GB

Additionally required: Timer / Clock (not included)

Der **Kombinations-Test** neuer Coronaviren (SARS-CoV-2-Antigene), Influenza A-Virus-Antigene, Influenza B-Virus-Antigene und RSV-Antigene

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
5s test

Product ID # **MF-71-5**

- 5 Test cassette
- 5 Sterile nasal swab
- 5 Prefilled sample Treatment Tube
- 1 Instruction manual D/GB

Additionally required: Timer / Clock (not included)

Der **Kombinations-Test** neuer Coronaviren (SARS-CoV-2-Antigene), Influenza A-Virus-Antigene, Influenza B-Virus-Antigene und RSV-Antigene

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25s test

Product ID # **MF-71-25**

- 25 Test cassette
- 25 Sterile nasal swab
- 25 Prefilled sample Treatment Tube
- 1 Workstation
- 1 Instruction manual D/GB

Additionally required: Timer / Clock (not included)

Der **Kombinations-Test** neuer Coronaviren (SARS-CoV-2-Antigene), Influenza A-Virus-Antigene, Influenza B-Virus-Antigene und RSV-Antigene

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
1. Disinfect your hands with a disinfectant after washing them.
2. Clean the surface on which the test will be performed.
3. Before testing, please read the operating instructions carefully.
4. Bring the test kit to room temperature (20-25°C) before use.
5. Twist off the lid of the extraction tube and remove the inner blue plug. (The purpose of the blue stopper is to prevent the product from leaking during transport. The blue plug must be removed before use!)
6. Insert the extraction tube into the opening* provided in the packaging to hold it in a vertical position.
7. Open the foil pouch of the test cassette, take out the test card and use it as soon as possible within 1 hour.

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Please note the perforation on the kits of the of the 1 and 5 self-test.

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
In the kit of the 25-series test for professional use a workstation is included.

Der **Kombinations-Test** neuer Coronaviren (SARS-CoV-2-Antigene), Influenza A-Virus-Antigene, Influenza B-Virus-Antigene und RSV-Antigene

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

Carefully remove the sterile nasal swab from the packaging. (Avoid touching the padded end).
Insert the swab into the left nostril to a depth of 2.5 cm from the edge of the nostril.

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Rotate the nasal swab 5 times on the nostril wall (mucosa) to ensure adequate sample collection.

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Repeat in the right nostril using the same nasal swab to obtain a sufficient sample from both nasal passages.

Der **Kombinations-Test** neuer Coronaviren (SARS-CoV-2-Antigene), Influenza A-Virus-Antigene, Influenza B-Virus-Antigene und RSV-Antigene

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

The **professional test** offers a choice of **three test variants**: Throat swab, nasopharyngeal swab or nasal swab.

- With the mouth wide open, wipe the base of the patient's tongue and pharyngeal tonsils with a swab. While doing so, gently rub back and forth, up and down on the posterior pharyngeal wall at least 3 times.

- Tilt the patient's head back and insert the swab through the nostril entrance and then slowly along the floor of the nasal cavity.
- When the swab tip reaches the posterior wall of the nasopharynx, gently rotate the swab a few times to collect as much secretion as possible. Remove carefully.

- Insert the swab into the nostril about 2 cm deep while carefully rotating it parallel to the palate until resistance is felt at the nostrils.
- Turn the swab several times against the nasal wall. Remove the swab, insert it into the other nostril and repeat the procedure.

Der **Kombinations-Test** neuer Coronaviren (SARS-CoV-2-Antigene), Influenza A-Virus-Antigene, Influenza B-Virus-Antigene und RSV-Antigene

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


Insert the swab with the swab sample into the extraction tube and break the swab at the predetermined breaking point. Leave the lower half in the tube and close the lid.

Squeeze the swab 10 times in the sample extraction tube. Then wait for 1 minute for the sample to react. Unscrew the lid on the top of the extraction tube to remove the blue plug underneath. Then screw the top cap back onto the extraction tube.

To perform the test, **two drops (approximately 60 µL)** of the sample solution are dropped into each sample well (S). There is the option of performing any single test or using all three sample variants.

Der **Kombinations-Test** neuer Coronaviren (SARS-CoV-2-Antigene), Influenza A-Virus-Antigene, Influenza B-Virus-Antigene und RSV-Antigene

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Before testing, carefully read the instructions for use and bring the test kit and components to a room temperature of 20 - 25°C.

The test procedure takes 15 minutes. Read the test results after 15 minutes.
If you read the test results after 20 minutes, the test results may be incorrect or invalid.
While you are waiting, do not touch the test card or pick it up from the tabletop.

When the test is complete, place the used components in a plastic bag. Seal it.
Dispose of it and disinfect your hands.

Der **Kombinations-Test** neuer Coronaviren (SARS-CoV-2-Antigene), Influenza A-Virus-Antigene, Influenza B-Virus-Antigene und RSV-Antigene

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
1. Positive for COVID-19 antigen or RSV: Two red lines, both the test line (T-line) and the control line (C line) indicate the red color.
NOTE: It does not matter if the line (T) is lighter or darker than the other; the result is "positive".
2. Positive for influenza A/B: Line A and the control line (C line) are colored to indicate positivity for influenza A. Line B and the control line (C line) are colored to represent positivity for influenza B. Line A, line B and the quality control line (C-line) indicate by their color that one is positive for both influenza A and influenza B.
NOTE: It does not matter if line A or line B is lighter or darker than the other two; the result is "positive".
3. Negative: a red stripe for the control line (C line).
4. invalid: the position of the control line (C line) in the observation window shows no color reproduction. This means that the test is invalid. You should repeat the sampling and test with new kits.

Der **Kombinations-Test** neuer Coronaviren (SARS-CoV-2-Antigene),
Influenza A-Virus-Antigene, Influenza B-Virus-Antigene und RSV-Antigene

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








The self-test is available in two packaging units:

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The professional test is available in available in one packaging unit:

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ZERTIFIKAT | CERTIFICATE | CERTIFICADO | CERTIFICATO | CERTIFICATE

Certificate
No. Q5 109172 0001 Rev. 00

Holder of Certificate: Shenzhen Microprofit Biotech Co., Ltd
Rm. 405, 406, Zone B 4/F
Rm. 205, 206-1, 207, West Side of Zone B/2/F
Haowei Building, No. 8 Langshan 2nd Road
Songgongshan, Songgongshan Community
Xili Street, Nanshan District
518057 Shenzhen
PEOPLE'S REPUBLIC OF CHINA

Certification Mark: 

Scope of Certificate: Design and Development, Production and Distribution of Immunochromatographic Assay Diagnostic Kit, Colloidal Gold Chromatographic Immunoassay Test Kit, and Dry-Type Immunofluorescence Quantitative Analyzer

The Certification Body of TÜV SÜD Product Service GmbH certifies that the company mentioned above has established and is maintaining a quality management system, which meets the requirements of the listed standard(s). All applicable requirements of the testing and certification regulation of TÜV SÜD Group have to be complied with. For details and certificate validity see: www.tuvsvud.com/ps-cert?cert=Q5_109172_0001_Rev_00

Report No.: GZ2043601

Valid from: 2021-03-24
Valid until: 2024-03-23

Date: 2021-03-24

Christoph Dicks
Head of Certification/Notified Body

Page 1 of 2
TUV SUD Product Service GmbH • Certification Body • Ridlenstraße 65 • 80339 Munich • Germany

ZERTIFIKAT | CERTIFICATE | CERTIFICADO | CERTIFICATO | CERTIFICATE

Certificate
No. Q5 109172 0001 Rev. 00

Applied Standard(s): EN ISO 13485:2016
Medical devices - Quality management systems - Requirements for regulatory purposes (ISO 13485:2016)
DIN EN ISO 13485:2016

Facility(ies): Shenzhen Microprofit Biotech Co., Ltd
Rm. 405, 406, Zone B 4/F, Rm. 205, 206-1, 207, West Side of Zone B/2/F, Haowei Building, No. 8 Langshan 2nd Road, Songgongshan, Songgongshan Community, Xili Street, Nanshan District, 518057 Shenzhen, PEOPLE'S REPUBLIC OF CHINA

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These documents can be requested at info@ultimed.de

Der **Kombinations-Test** neuer Coronaviren (SARS-CoV-2-Antigene), Influenza A-Virus-Antigene, Influenza B-Virus-Antigene und RSV-Antigene

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 TO-04-005 A/2

Declaration of Conformity

Manufacturer: **Shenzhen Microprofit Biotech Co., Ltd.**
Rm. 405, 406, Zone B 4/F, Rm. 205, 206-1, 207, West Side of Zone B/ 2F, Haowei Building, No. 8 Langshan 2nd Road, Songpingshan, Songpingshan Community, Xili Street, Nanshan District, Shenzhen, P.R. China

European Representative: **CMC MEDICAL DEVICES & DRUGS, S.L.**
C/ Horacio Lengo n18 - C.P 29006 - Málaga-Spain

Product Name: **fluorecare SARS-CoV-2 & Influenza A/B & RSV Antigen Combo Test Kit**

Common Name: **SARS-CoV-2 & Influenza A/B & RSV Antigen Combo Test Kit (Colobal Gold Chromatographic Immunoassay)**

Brand: **fluorecare®**

Catalogue No.: **MF-71-1, MF-71-2, MF-71-5**

Classification: **Self-testing Device of IVDD 98/79/EC**

Conformity Assessment Route: **Annex III of IVDD 98/79/EC**

| | | |
|--------------------------|---|--|
| STANDARDS APPLIED | EN 13612:2002/AC: 2002 EN ISO 14971:2012 EN ISO 18113-1:2011 EN ISO 15223-1:2016 | EN ISO 13485:2016 EN ISO 23640:2015 EN ISO 18113-2:2011 EN 13641:2002 |
|--------------------------|---|--|

We the manufacturer herewith declare on our solo responsibility that the above mentioned products meet the provisions of the following EC Council Directives and Standards. All supporting documentations are retained under the premises of the manufacturer.
The products comply with the essential requirements in accordance with Annex I of the IVDD 98/79/EC.

DIRECTIVES

General applicable directives:
In Vitro Diagnostic Medical Device Directive: **COUNCIL DIRECTIVE 98/79/EC of 27 October 1998 concerning in vitro diagnostic medical devices (IVDD 98/79/EC).**

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CeCert.

CERTIFICATE

DIRECTIVE 98/79/EC
EC DESIGN-EXAMINATION

CeCert Sp. z o.o. hereby confirms that manufactured by

Shenzhen Microprofit Biotech Co., Ltd.
Rm. 405, 406, Zone B/4F, Rm. 205, 206-1, 207, West Side of Zone B/ 2F, Haowei Building, No. 8 Langshan 2nd Road, Songpingshan, Songpingshan Community, Xili Street, Nanshan District, Shenzhen, P.R. China

in vitro diagnostic medical device for self-testing

fluorecare SARS-CoV-2 & Influenza A/B & RSV Antigen Combo Test Kit
catalogue numbers: **MF-71-1, MF-71-2, MF-71-5**

in term of the design conforms to the requirements of Annex III section 6 to Directive 98/79/EC (as amended) implemented into Polish Law, as evidenced by the assessment conducted by CeCert Sp. z o.o.

CE 2934

Validity date: 12.05.2022 - 26.05.2025
Edition issue date: 18.05.2022

Check it 

Kamil Szczurowski
Director of *in Vitro* Diagnostic Medical Device Certification Department

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DECLARATION OF CONFORMITY

MANUFACTURER Shenzhen Microprofit Biotech Co., Ltd.
Rm. 405, 406, Zone B /4F, Rm. 205, 206-1, 207, West Side of Zone B/ 2F, Haowei Building, No. 8 Langshan 2nd Road, Songpingshan, Songpingshan Community, Xili Street, Nanshan District, Shenzhen, P.R. China

EUROPEAN REPRESENTATIVE CMC MEDICAL DEVICES & DRUGS, S.L.
C/ Horacio Lengo n18 - CP 29006 - Malaga-Spain

PRODUCT See the attachment

CLASSIFICATION Other Device of IVDD 98/79/EC

CONFORMITY ASSESSMENT Annex III of IVDD 98/79/EC

We the manufacturer herewith declare on our solo responsibility that the above mentioned products meet the provisions of the following EC Council Directives and Standards. All supporting documentations are retained under the premises of the manufacturer. The products comply with the essential requirements in accordance with Annex I of the IVDD 98/79/EC.

| | | |
|--------------------------|--|--|
| STANDARDS APPLIED | EN 13612:2002/AC: 2002 EN ISO 14971:2012 EN ISO 18113-1:2011 EN ISO 15223-1:2016 EN 62366-1:2015 | EN ISO 13485:2016 EN ISO 23840:2015 EN ISO 18113-2:2011 EN 13641:2002 |
|--------------------------|--|--|

PLACE Shenzhen, China

DATE OF ISSUE 2022-4-23

SIGNATURE   
General Manager

ATTACHMENT


| Catalog No. | Product name |
|-------------|--|
| MF-59 | SARS-CoV-2 Antibody Test Kit (Colloidal Gold Chromatographic Immunoassay) |
| MF-60 | SARS-CoV-2 Spike Protein Test Kit (Colloidal Gold Chromatographic Immunoassay) |
| MF-61 | SARS-CoV-2 IgG and IgM Antibody Combined Test Kit (Colloidal Gold Chromatographic Immunoassay) |
| MF-62 | SARS-CoV-2 Antibody Test Kit (Fluorescence Immunoassay) |
| MF-63 | SARS-CoV-2 Spike Protein Test Kit (Fluorescence Immunoassay) |
| MF-64 | SARS-CoV-2 IgG and IgM Antibody Combined Test Kit (Fluorescence Immunoassay) |
| MF-66 | fluorecare COVID-19 & Influenza A/B Antigen Combo Test Kit(Colloidal Gold) |
| MF-67 | SARS-CoV-2 Antigen Test Kit (Fluorescence Immunoassay) |
| MF-68 | SARS-CoV-2 Antigen Test Kit (Colloidal Gold Chromatographic Immunoassay) |
| MF-69 | SARS-CoV-2 & Influenza A/B Antigen Combo Test Kit (Colloidal Gold Chromatographic Immunoassay) |
| MF-70 | SARS-CoV-2 Neutralizing Antibody Test Kit (Fluorescence Immunoassay) |
| MF-71 | SARS-CoV-2 & Influenza A/B & RSV Antigen Combo Test Kit (Colloidal Gold Chromatographic Immunoassay) |
| MF-90 | SARS-CoV-2 Neutralizing Antibody Test Kit (Colloidal Gold Chromatographic Immunoassay) |
| MF-91 | SARS-CoV-2 Antigen Test Kit (Colloidal Gold Chromatographic Immunoassay) (Saliva) |

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 These document can be requested at info@ultimed.de

Der **Kombinations-Test** neuer Coronaviren (SARS-CoV-2-Antigene),
Influenza A-Virus-Antigene, Influenza B-Virus-Antigene und RSV-Antigene

Der Kontakt | The contact

MICROPROFIT

Shenzhen Microprofit Biotech Co., Ltd.

Songpingshan Community, Xili Street, Nanshan District,
Shenzhen, P.R. China

CMC MEDICAL DEVICES & DRUGS, S.L.
C / Horacio Lengo n 18 · C.P 29006 · Málaga-Spain

Importeur **ulti med Products (Deutschland) GmbH**

Reeshoop 1,

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www.ultimed.org

E-Mail: info@ultimed.org / sales@ultimed.org