SIMPLE/STICK Flu A+B

USO EXCLUSIVO PROFESIONAL

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SIMPLE/STICK Flu A+B

FOR PROFESSIONAL USE ONLY

PLEASE READ THESE INSTRUCTIONS CAREFULLY BEFORE USING THE TEST.

INTRODUCTION

• Intended use

A one-step immunochromatographic test to simultaneously qualitatively detect Influenza A, and Influenza B antigens in human respiratory samples (nasopharyngeal and oropharyngeal swabs, nasopharyngeal washes and aspirates) on separate strips.

• General information

Influenza virus tends to appear during the winter months, with incidence highest for around a 5 month period, before it completely disappears when spring arrives¹.

The Flu A+B test can detect Influenza A and B virus nucleoproteins in the same sample.

This product is exclusively for *in vitro* diagnostic use.

The results of this test are intended to support data derived from a clinical evaluation of the patient; the test does not constitute a definitive diagnosis in itself. In the event of obtaining a negative result for a patient presenting the typical symptomatology from these types of infections, confirming the result using RT-PCR is recommended.

This is a test for professional use. This is not a self-diagnostic test.

Intended test population

The test is intended for use on the general public, as any individual is susceptible to becoming infected with the virus.

Disease or infection incidence among the target population for this product

Worldwide, from the start of November 2019 to the end of December 2020, 614,907 respiratory samples collected during hospital and outpatient surveillance were positive for Influenza (19% of all analysed samples). Of these positive samples, 63% and 37% were typed as influenza A and influenza B respectively.² The WHO estimates that annual flu epidemics cause approximately 3 to 5 million cases of serious disease, and 300,000 to 500,000 deaths.

• Characteristics of the virus and its infection

Influenza A viruses are divided into subtypes based on two proteins on the virus' surface: hemagglutinin (H) and neuraminidase (N). The influenza A virus subtypes

currently circulating more commonly among people are: A (H1N1) and A (H3N2). Influenza B viruses, on the other hand, are classified into two lineages: B/Yamagata and B/Victoria.

The Influenza virus is an RNA virus that presents the hemagglutinin (HA) and neuraminidase (NA) proteins in its envelope. The HA protein is responsible for the interaction with host cells, and therefore, the virus' entry, while the NA protein participates in both virus replication and its exit from the cell. The HA and NA viral proteins show most variation in terms of antigens, while the virus' nucleoprotein, detected by the antibodies used for both Influenza A and Influenza B, is much more conserved.

Symptoms associated with flu virus infection vary from a mild respiratory disease, limited to the upper respiratory tract - characterised by fever, throat pain, nasal secretion, cough, headache, muscle pain and fatigue - to severe pneumonia, which in some cases can prove fatal due to the Influenza virus itself or a secondary bacterial infection in the lower respiratory tract. In some cases, infection by Influenza virus can also cause a wide range of non-respiratory complications that affect the heart, central nervous system, and other organ systems.³

THE BASIS OR BASIC PRINCIPLES OF THE TEST

The Flu A+B test is an immunochromatographic test with one strip inside a housing. The strip uses a combination of:

a. Purple latex particles conjugated to an antigen recognised by a specific antibody for said antigen bound to the membrane, comprising the so-called test control line (position C on the housing).

b. Blue latex particles conjugated to a specific antibody against the Influenza type B virus nucleoprotein that cooperates with another specific antibody for said protein located on the membrane, below the control line (position T2 on the housing).

c. Red latex particles conjugated to a specific antibody against the Influenza type A virus nucleoprotein that cooperates with another specific antibody for said protein located on the membrane, below the line that detects Influenza type B (position T1 on the housing).

This test first involves treating the sample with the sample dilution buffer to extract the Influenza virus (included in this kit). After extraction, only a certain volume of the extract needs to be added to the reactive strips, before waiting for 15 minutes.

When the extracted sample flows across the strip, the coloured particles migrate. In the event of a positive sample, the antibodies present in the particles will react with the virus antigens, forming a particle-antigen complex, which will also be retained on the membrane by other specific antibodies for that antigen

Different coloured lines will appear depending on the virus content in the sample. These lines are used to interpret the result after 15 minutes of incubation at room temperature (see Fig. 1).

These instructions for use apply for all commercial product references: 9.080.XXX.YY.ZZZ.

MATERIAL INCLUDED IN THE KIT

Simple Format

- 20 reaction devices.
- 1 vial with a dropper containing the sample dilution buffer.
- Instructions for use.

Depending on the reference, the kit could include also:

- 20 sterile swabs
FLOQ Swabs Copan ref. 502CS01  0123
- 20 disposable plastic pipettes.
- 20 x 1.5 mL graduated plastic microtubes with a lid.
- 1 vial with a dropper containing the sample dilution buffer.
- Positive control and corresponding instructions for use.
- Datamatrix

Stick Format

- 20 reaction devices.
- 1 vial containing the sample dilution buffer.
- 20 test tubes.
- 3 stands for test tubes.
- Instructions for use.

The product's safety sheet contains information about the composition of the reagents. An electronic copy of the sheet can be requested via email to msds@operon.es

MATERIAL NOT INCLUDED IN THE KIT

- Chronometer

Depending on the reference the kit could not include:

- Micropipette
- Plastic microtube.

Controls are available for the user to validate the results, which can also be obtained as an independent commercial reference.

EQUIPMENT FOR INTERPRETING THE RESULTS

The test can be interpreted visually or using OPERON's immunochromatographic test reader.

For this diagnosis, there is evidence of a 99% agreement between visual interpretation and the OPERON IC strip reader.

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PRECAUTIONS

Check all the applicable precautions from this list example:

1. Patient samples must be handled with care as they may contain infectious agents. Disposable gloves and other protective measures that may eventually be required must be used throughout the process.
2. The sample dilution buffer contains sodium azide as an anti-microbial agent. Avoid direct contact with skin and mucous membranes. Dispose of appropriately. Do not use the buffer if there are any signs of contamination or precipitation.
3. Do not store or prepare food, eat, drink, or smoke in the area where the reagents and samples are handled.
4. On completing the task, clean all work surfaces with soap and water and disinfect them with a suitable solution. Lastly, dispose of the protective gloves correctly and wash your hands with soap and water, rubbing them well.
5. Do not swap the components of kits with different batch numbers.
6. If the test is stored refrigerated, allow all the kit's components and the samples to reach room temperature, as cold reagents and/or samples can reduce the test's functionality. Waiting for 20 to 30 minutes to reach room temperature is recommended.
7. Only use the reagents *in vitro*.
8. The immunochromatographic strips are single use. A strip is used per sample to be evaluated.
9. Do not use any of the kit's components after the use-by-dates.
10. If the outer box is damaged, the product can still be used providing none of the components have been damaged.
11. Throw the test away if the primary packaging is either broken or impaired.
12. Remove the strip or device from the aluminium bag when the test is going to be performed to avoid unnecessary excessive exposure of the product to any potentially damaging environmental factors.
13. Adding the correct extracted sample volume to sample addition zone indicated on the housing is very important. If lower than indicated, the chromatography may not process correctly because insufficient sample reaches the reaction zone. If higher than indicated, the chromatography may not process correctly due to the sample volume surpassing the strip's absorption capacity.
14. Please dispose of the used product in accordance with current legislation.
15. Do not use the test if there are any coloured lines in the results zone prior to use.
16. Mix all respiratory samples thoroughly before performing the test in order to ensure a representative sample is used.
17. When wash/aspirate type samples are analysed, avoid sampling from this zone if the sample still has viscous zones after mixing. If the pipette blocks, return the sample to the container and sample it again.
18. Do not use sputum or saliva samples as invalid results may be obtained.

19. Do not discard the kit's outer box until all the content has been used. The box contains essential information about the product's CE marking and the batch number.

20. Any serious incident involving the product must be reported to the manufacturer and the relevant authority in the member state where the user and/or patient are located.

STORAGE

The Simple FluA+B product can be stored at any temperature from +2 to +30 °C.

The use-by date of each component is printed on the aluminium wrappers.

The reactive strips must be used once at room temperature (if stored refrigerated) after the protective packaging has been opened.

The product maintains the claimed stability once open.

SAMPLES

- The test is designed to analyse human respiratory samples from nasopharyngeal and oropharyngeal swabs, nasopharyngeal washes and aspirates.
- Samples resuspended in the product's sample dilution buffer must be analysed immediately after collection, it is not recommended to store them for later use. For samples collected in a viral transport medium different to the product's sample dilution buffer, please see the specifications for each medium when it comes to the stability of resuspended samples and storage conditions.
- Limit sample freezing/thawing cycles, given they could negatively affect the results.
- Do not centrifuge the samples prior to their use with the Simple Flu A+B test, given that the elimination of cell material could negatively affect the test's sensitivity.
- There is evidence that the nasopharyngeal washes and aspirates are more effective for detecting viruses than nasopharyngeal/oropharyngeal swab samples.⁴
- Nasopharyngeal wash/aspirate type samples must be mixed very well prior to taking the aliquot for analysis. If using dry swabs, ensure that they resuspend correctly in the dilution buffer, as indicated later.
- Several swab types have been evaluated for sample collection. If it is a oropharyngeal sample, the use of oropharyngeal Flocked Swab from Copan Flock Technologies included in the kit is recommended. For nasopharyngeal samples, the use of Goodwood nasal/nasopharyngeal swabs (ref. GW-1237NP) or Copan Flock Technologies FLOQSwabs (ref. 503CS01) is recommended.
- Samples have been analysed in different transport mediums. The test is compatible with at least the following: phosphate-buffered saline (PBS), viral transport medium from Vircell and CPM Scientifica. If using different transport mediums to the product's dilution buffer, the test's characteristics may vary depending on the type and volume used.

PREPARING SAMPLES AND THE TEST PROCEDURE

General note: disposable gloves must be used throughout the test due to the handling of infectious samples. On completion, follow the hygiene procedures in point 4 of the "Precautions" section.

It is very important to mix nasopharyngeal washes and aspirate samples well before taking the aliquot for analysis (a vortex mixer can be used, or failing that, stir vigorously by hand).

Mix oropharyngeal/nasopharyngeal swabs well to ensure swab sample extraction (if present) and to ensure the aliquot for analysis is mixed.

Do not centrifuge under any circumstances, as eliminating cell material could adversely affect the test's sensitivity.

Allow any refrigerated or frozen samples to reach room temperature before analysing them.

THE PROCEDURE WILL VARY DEPENDING ON BOTH THE SPECIMEN TO BE ANALYSED AND THE FORMAT OF THE TEST (Simple or Stick).

▪ **Procedure for freshly collected oropharyngeal/nasopharyngeal swabs without a transport medium.**

SIMPLE FORMAT

1. Dispense the sample dilution buffer into the plastic microtube supplied with the kit up to the third mark from the bottom of the tube (between the 0.5 mL and 1.0 mL marks), corresponding to 0.75 mL.

In the case of using a commercial reference that does not include this accessory, dispense 0.75 mL with a micropipette (not included in the kit) in a 1.5 mL plastic microtube.

2. If it is a oropharyngeal sample, collect the sample using the swabs included in the kit. For nasopharyngeal samples, using swabs reported in "SAMPLE" section is recommended or swabs with similar characteristics, should none of these be available.

3. Introduce the swab soaked in sample into the microtube containing the dilution buffer.

4. Vigorously twist the swab in the buffer for 30 seconds to extract the sample correctly.

5. Remove the swab, pressing it against the microtube walls to ensure as much sample as possible is released.

6. Dispose of the swab.

7. Remove the reaction device from its packaging and place it on a flat surface.

8. Add **7 drops** of the extracted sample to the sample addition zone of the strip (rectangular window indicated with an arrow) using the plastic pipette supplied with the kit. Add slowly drop by drop to ensure the sample absorbs correctly.

In the case of using a commercial reference that does not include this accessory, add 0.14 mL to the test with a micropipette (not included in the kit). Pay attention in adding slowly drop by drop to ensure the sample absorbs correctly.

9. Wait **15 minutes**, then read and interpret the result.

▪ Procedure for:

→ **oropharyngeal/nasopharyngeal swab samples in viral transport medium**

→ **nasopharyngeal washes and aspirates.**

SIMPLE FORMAT

1. Dispense 150 μL of sample dilution buffer into a microtube.

2. Add 150 μL of sample using a pipette and gently mix, ensuring to avoid any foam being generated (avoid using a vortex mixer).

3. Remove the reaction device from its packaging and place it on a flat surface.

4. Add 140 μL of the mix (sample+dilution buffer) to the sample addition zone of the strip (rectangular window indicated with an arrow). Add slowly, drop by drop, to ensure the sample absorbs correctly.

5. Wait **15 minutes**, then read and interpret the result.

STICK FORMAT

1. Dispense 150 μL of sample dilution buffer into a test tube.

2. Add 150 μL of sample using a pipette and gently mix, ensuring to avoid any foam being generated (avoid using a vortex mixer).

3. Remove the reaction device from its packaging and place it inside the test tube in a vertical position.

4. Wait **15 minutes**, then read and interpret the result.

READING THE RESULTS

The strips shown in Fig. 1 are an example of the different results that can be obtained with the Flu A+B test.

Different coloured lines can appear on each reagent strip, bordered by black horizontal lines printed on one side of the housing:


a. Purple line: this is the control line, which should always appear, as it indicates that the chromatography worked correctly.

b. Blue line: this indicates that the sample contains influenza type B virus.

c. Red line: this indicates that the sample contains influenza type A virus.

The purple control line must always appear.

Figure 1 shows examples of possible results.

Strip 1: **NEGATIVE** result

Only a **PURPLE** horizontal line appears on both strips aligned with the letter “C” on the housing. This is the control line, which should always appear, as it indicates that the chromatography worked correctly.

Strips 2-4: **POSITIVE** results

Strip 2: Influenza A Positive + Influenza B Negative

A **PURPLE** control line appears, aligned with the letter “C” and a **RED** line appears aligned with the letter “T1”.

Strip 3: Influenza A Negative + Influenza B Positive

A **PURPLE** (control) line appears aligned with the letter “C” and a **BLUE** line appears aligned with the letter “T2”.

Strip 4: Influenza A Positive + Influenza B Positive

A **PURPLE** (control) line appears aligned with the letter “C”, a **RED** line appears aligned with the letter “T1” and a **BLUE** line appears aligned with the letter “T2”.

Strip 5: **INVALID** results

The control lines do not appear or the purple colour of the control lines is clearly altered. This is an indication that the test has not worked correctly. In the event of an invalid result, repeating the test using a new reaction device is recommended, strictly adhering to the instructions for use described in this manual.

Any line that appears after 15 minutes of reaction time will have no diagnostic value.

NOTE: a doctor is responsible for the final and definitive diagnosis. This test only detects respiratory influenza viruses in a sample, it cannot be used to definitively conclude that the person is infected with those viruses.

The test can be interpreted visually or using OPERON’s immunochromatographic test reader. The quick guide for interpreting the test using that reader is available from the following website address: <https://operon.es/operon-lateral-flow-reader/>

Please note: check that no dust particles, fibres, etc. are present in the results window prior to introducing the test into the reader, as they can interfere with the reading.

QUALITY CONTROL

The test is invalid should no (coloured) lines appear, either because it was performed incorrectly or because the reagents have deteriorated. In this case, repeat the analysis, strictly following the operating protocol detailed in these instructions for use.

WARNING: Including our controls with a known result is recommended to ensure that correct data is obtained.

LIMITATIONS OF THE PROCEDURE

1. The Flu A+B test can (separately) identify Influenza A and B, detecting their presence in human respiratory samples such as nasopharyngeal/oropharyngeal swabs and nasopharyngeal wash/aspirates, as long as the viral load is equal to or above the assay’s limit of detection.

2. This test is qualitative, not quantitative, although the intensity of the positive bands does relate to the amount of detectable analyte in the sample.

3. A large amount of respiratory samples were evaluated to ensure the test functions correctly. The results correlated well against other reference techniques (such as RT-PCR Real Time). However, this study does not rule out any potential interference to the test’s performance on analysing other respiratory samples.

4. Monitoring the reaction time is important. If the reaction time is shorter than indicated, samples containing a quantity of analyte above the sensitivity limit will be clear on observation, but those on the limit will not appear. If the reaction time is longer than indicated, the test’s sensitivity may be impaired and give rise to incorrect interpretations.

5. The test results should be interpreted alongside the available data from epidemiological studies, the patient’s clinical evaluation, and other diagnostic procedures.

6. Adding the correct amount of sample to the reactive strip (7 drops with the plastic pipette supplied in the kit or 140 μL) is very important, as a smaller amount may cause false negative results, whereas a larger amount could prevent the chromatogram from developing correctly.

7. Unsuitable collection, transportation, and/or handling of the sample could impair the results.

8. In the case of samples with a very high Influenza B viral load, a weaker, blue, non-specific signal could appear in correspondence with the Flu A line. This does not mean that there is infection by the Influenza A virus.

9. The antibodies present in the test may not detect all antigenic variants of the influenza viruses, or any new

strains that appear over time. The Analytical Sensitivity section lists all the strains analysed with the Flu A+B test.

ANALYTICAL SENSITIVITY

▪ Type A

- Analytical sensitivity with Influenza type A (H1N1), strain A/New Caledonia/20/99 IVR116 chemically inactivated: 8 ng/mL.

- Analytical sensitivity with Influenza type A (H1N1), strain A/New Caledonia/20/99 heat inactivated: 19.5 TCID₅₀/mL.

- Analytical sensitivity with Influenza type A (H1N1), strain A/Wisconsin/67/2005 (viable virus): 3,6 x 10⁵ CEID₅₀/mL.

▪ Type B

- Analytical sensitivity with Influenza type B, lineage Yamagata, strain B/Florida/07/04 chemically inactivated: 16 ng/mL.

- Analytical sensitivity with Influenza type B, lineage Yamagata, strain B/Florida/07/04 heat inactivated: 2.5 x 10³ TCID₅₀/mL.

- Analytical sensitivity with Influenza type B, lineage Yamagata, strain B/Florida/07/04 (viable virus): 3,6 x 10⁵ CEID₅₀/mL.

The Influenza strip analytical sensitivity was analysed with other variants of the virus (A and B):

Strain	Subtype/ lineage	Analytical sensitivity
A/New Caledonia/20/99	H1N1	8 ng/mL
A/Solomon Island/03/06	H1N1	31 ng/mL
A/Beijing/262/95	H1N1	31 ng/mL
A/Taiwan/1/86	H1N1	16 ng/mL
A/Brisbane/10/07	H3N2	8 ng/mL
A/Kiev/301/94	H3N2	16 ng/mL
A/Panama/2007/99	H3N2	16 ng/mL
A/Shangdong/9/93	H3N2	16 ng/mL
A/Wisconsin/67/05	H3N2	31 ng/mL
A/Florida/04/06	Yamagata	8 ng/mL
B/Hong Kong/5/72	Victoria	31 ng/mL
B/Tokyo/53/99	Victoria	16 ng/mL
B/Victoria/504/00	Victoria	31 ng/mL

ANALYTICAL SPECIFICITY

The monoclonal antibodies used to detect Influenza A and B specifically detect a highly conserved epitope from the nucleoprotein of each virus.

DIAGNOSTIC SENSITIVITY AND SPECIFICITY

The Flu A+B test was evaluated by retrospectively analysing different respiratory sample specimens.

The reference techniques were:

- RT-PCR CLART® PneumoVir 2 (Genomica).
- PCR Cobas Liat Influenza A/B (Roche).
- Aries® Real Time PCR Flu A/B&RSV (Luminex).
- Real SARS-CoV-2/Flu/RSV (Operon).

The results are shown below and separated by analyte, specimen, and reference technique. 190 respiratory samples were analysed for Influenza A and 77 of these samples were clinically diagnosed as positive. 145 respiratory samples were analysed for Influenza B and 14 of these samples were clinically diagnosed as positive.

Oropharyngeal swab samples in viral transport medium characterised by RT-PCR CLART® PneumoVir 2 (Genomica)

36 samples negative for Influenza A/B were analysed, along with 34 positive for Influenza A and 6 positive for Influenza B.

Influenza A		RT-PCR CLART® PneumoVir 2	
		Positive	Negative
Flu A+B test	Positive	20	0
	Negative	14	36
Sensitivity agreement = 58,8% Specificity agreement ≥ 99,9% PPV: 100 % NPV: 72,0%			

Influenza B		RT-PCR CLART® PneumoVir 2	
		Positive	Negative
Flu A+B test	Positive	2	0
	Negative	2	36
Sensitivity agreement = 50,0% Specificity agreement ≥ 99,9% PPV: 100% NPV: 94.7%			

Oropharyngeal swab samples in viral transport medium characterised by RT-PCR Real Time Cobas Liat Influenza A/B (Roche)

26 samples negative for Influenza A/B and 22 positive for Influenza A were analysed.

Influenza A		RT-PCR Cobas Liat Influenza A/B	
		Positive	Negative
Flu A+B test	Positive	14	1
	Negative	8	25
Sensitivity agreement = 63,6% Specificity agreement = 96,1% PPV: 93,6% NPV: 75,8%			

Influenza B		RT-PCR Cobas Liat Influenza A/B	
		Positive	Negative
Flu A+B test	Positive		0
	Negative		48
Specificity agreement ≥ 99.9%			

Oropharyngeal swab samples characterised by Real SARS-CoV-2/Flu/RSV (Operon)

47 samples negative for Influenza A/B were analysed.

Influenza A		RT-PCR	
		Positive	Negative
RSV-Flu A+B test	Positive		1
	Negative		46
Specificity agreement = 97.9% PPV: N/A NPV: 100%			

Influenza B		RT-PCR	
		Positive	Negative
RSV-Flu A+B test	Positive		0
	Negative		47
Specificity agreement ≥ 99.9% PPV: N/A NPV: 100%			

Nasopharyngeal aspirate samples characterised by Luminex Aries® Real Time PCR Flu A/B&RSV

Influenza A: 73 samples were analysed, 21 of them originally characterized as positive.

Influenza B: 71 samples were analysed, 10 of them originally characterized as positive.

Influenza A		RT-PCR	
		Positive	Negative
Flu A+B test	Positive	18	1
	Negative	3	51
Sensitivity agreement = 85,7% Specificity agreement = 98,1% PPV: 94,7% NPV: 94,4%			

Influenza B		RT-PCR	
		Positive	Negative
Flu A+B test	Positive	9	0
	Negative	1	61
Sensitivity agreement = 90,0% Specificity agreement ≥ 99.9% PPV: 100% NPV: 98,4%			

REPEATABILITY

Five replicates of 1/2 serial dilutions of characterised internal standards were analysed (sensitivity curve), along with real samples established as NC ("negative sample"), LPC ("weak positive sample") and PC ("positive sample") for each analyte detected by the test, measured on the same day by the same person. High repeatability was obtained with differences between the replicates of the standard for each analyte equal to or below one 1/2 dilution and the results with the real samples were the same.

REPRODUCIBILITY

• INTER-DAY PRECISION

Over five consecutive days and in duplicate, the same Simple Flu A+B test batch was measured for a sensitivity curve and a group of real samples established as NC ("negative sample"), LPC ("weak positive sample") and PC ("positive sample") for each of the analytes detected using the test. Sensitivity differences equal or less than one 1/2 dilution were obtained and the results with the real samples were the same.

• INTER-OPERATOR PRECISION

In triplicate, three persons measured a sensitivity curve a group of real samples established as NC ("negative sample"), LPC ("weak positive sample") and PC ("positiva sample") for each one of the analytes detected with the test. Sensitivity differences equal or less than one 1/2 dilution were obtained and the results with the real samples were the same.

• INTER-BATCH PRECISION

In parallel and in duplicate using three different test batches, a sensitivity curve was evaluated and a group of real samples established as NC ("negative sample"), LPC ("weak positive sample") and PC ("positive sample") for each of the analytes detected with the test. The analysis was performed by one person on the same day. Sensitivity differences equal or less than one 1/2 dilution were obtained and the results with the real samples were the same.

THE HOOK OR PROZONE EFFECT

It was confirmed that the Flu A+B test:

1. Does not present a Hook effect with Influenza A/Virginia/ATCC1/2009 (Ref. VR-1736 from ATCC) until a viral infection load of 3.2×10^5 TCID₅₀/mL.

2. Does not present a Hook effect with Influenza B/Florida/04/06 (ref. VR-1804 de ATCC) until a viral infection load of $2,2 \times 10^7$ CEID₅₀/mL.

INTERFERING SUBSTANCES

The following substances showed no effect on the test results when added at the indicated concentrations:

Active ingredient	Concentration
Blood	0.5 %
Mucin	0.5 mg/mL
Benzocaine	1 mg/mL
Hyaluronic acid and allantoin	1 % v/v
Phenylephrine	15 % v/v
Oxymetazoline	15 % v/v
Mometasone	0.01 mg/mL
Fluticasone propionate	2.5 % v/v
Dichlorobenzyl alcohol and amylmetacresol	5% w/v (50 mg/mL)
Sodium cromoglycate	15 % v/v
Oseltamivir	5 mg/mL
Mupirocin	10 µg/mL
Tobramycin	4 µg/mL
Amoxicillin	1 mg/mL
Biotin	1 µg/mL
Alkalol	1:10
Phenolic compounds	5 % v/v
Ibuprofen	2.2 mg/mL
Paracetamol	1 mg/mL
Zanamivir	5 mg/mL
Ribavirin	1 mg/mL
Acetylsalicylic acid	0.3 mg/mL
Budesonide	1 % v/v

CROSS-REACTIVITY

Evidence for no cross-reactivity with the following micro-organisms has been obtained:

Micro-organism	Analysed concentration
Bordetella pertussis	1×10^6 CFU/mL
Haemophilus influenzae	1×10^6 CFU/mL
Legionella pneumophila	1×10^6 CFU/mL
Inactivated mycobacterium tuberculosis	1:10 Dilution (0.27 g/mL)
Streptococcus pneumoniae	1×10^6 CFU/mL
Streptococcus salivarius	1×10^6 CFU/mL
Streptococcus pyogenes	1×10^6 CFU/mL
Staphylococcus aureus	1×10^6 CFU/mL
Staphylococcus epidermis	1×10^6 CFU/mL
Candida albicans	1×10^6 CFU/mL

Pseudomonas aeruginosa	1×10^6 CFU/mL
Chlamydia pneumoniae	7.5×10^6 IFU/mL
Mycoplasma pneumoniae	1.43×10^6 CFU/mL
Nasal wash	N/A
Inactivated coronavirus 229E	0.94×10^4 TCID ₅₀ /mL
Inactivated coronavirus OC43	5×10^4 TCID ₅₀ /mL
Viable coronavirus OC43	0.89×10^4 TCID ₅₀ /mL
Inactivated coronavirus NL63	1.7×10^4 TCID ₅₀ /mL
Inactivated MERS	1.7×10^4 TCID ₅₀ /mL
Inactivated adenovirus 7A	1.5×10^5 TCID ₅₀ /mL
Metapneumovirus TN/91-316	2.8×10^5 TCID ₅₀ /mL
Parainfluenza 1 HPIV1/FRA/27344044/2007	0.89×10^4 TCID ₅₀ /mL
Parainfluenza 2 (Greer)	0.89×10^5 TCID ₅₀ /mL
Parainfluenza 3 (NIH 47885)	1.6×10^6 TCID ₅₀ /mL
Parainfluenza 4B (19503)	5×10^5 TCID ₅₀ /mL
Rhinovirus 48 (1505)	5×10^5 TCID ₅₀ /mL
Enterovirus D 68 (USA/2018-23089)	2.8×10^5 TCID ₅₀ /mL
Inactivated SARS-CoV-2	6.5×10^5 TCID ₅₀ /mL
RSV B1	$8,9 \times 10^4$ TCID ₅₀ /mL

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OPERON, S.A.. - Camino del Plano, 19 -E-50410 Cuarte