



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	173CS05
GENERAL DESCRIPTION	SPREADERS A T B/05PCS
DETAILED DESCRIPTION	1 box containing 200 paper-film pouches each containing 5 T-shape spreaders

PRODUCT COMPONENTS

SPREADERS	SPREDERS A T	DIMENSIONAL SPECIFICATIONS
		TOTAL LENGTH: 141.5 mm \pm 0.5 mm DIAMETER: 4.0 mm \pm 0.2 mm THICKNESS "T": 4.0 mm \pm 0.2 mm LENGTH "T": 35.0 mm \pm 0.5 mm

PACKAGING SPECIFICATION

PACKAGING	PAPER-FILM POUCH
BOX QTY	1000
ENVELOPE QTY	5

IMPORTANT INFO

STERILIZATION METHOD	GAMMA IRRADIATION
STORAGE TEMPERATURE	2/30°C – DO NOT FREEZE
SHELF LIFE	36 MONTHS