



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	176CS20
GENERAL DESCRIPTION	NEEDLES B/20PCS 4000PKG
DETAILED DESCRIPTION	1 outer box contain 4 inner boxes each contains 50 paper-film pouches each contains 20 needles

PRODUCT COMPONENTS

		SPECIFICATIONS
SPREADERS	NEEDLE	TOTAL LENGTH: 199.00 mm ± 1.00 mm TIP DIAMETER: 1.50 mm ± 0.20 mm MATERIAL: POLYSTYRENE

PACKAGING SPECIFICATION

INNER BOX QTY	1000
OUTER BOX QTY	4000
PACKAGING	PAPER-FILM POUCH
ENVELOPE QTY	20

IMPORTANT INFO

STERILIZATION METHOD	GAMMA IRRADIATION
STORAGE TEMPERATURE	2/30°C – DO NOT FREEZE
SHELF LIFE	36 MONTHS

This bulletin (Rev.01-30/06/2020) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.