



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	178CS10
GENERAL DESCRIPTION	ANSE PER INOC 1UL SOFT 2000PKG
DETAILED DESCRIPTION	1 outer box contain 4 inner boxes each composed of 50 paper-film pouches each containing 10 1uL soft loops.

PRODUCT COMPONENTS

		SPECIFICATIONS
LOOP	ANSA 1uL SOFT LIGHT GREEN	TOTAL LENGTH: 199.00 mm \pm 1.00 mm EXTERNAL DIAMETER RING: 3.00 mm \pm 0.20 mm INTERNAL DIAMETER RING: 0.80 mm \pm 0.20 mm THICKNESS RING: 1.10 mm \pm 0.20 mm MATERIAL: POLYSTYRENE

PACKAGING SPECIFICATION

INNER BOX QTY	500
OUTER BOX QTY	2000
PACKAGING	PAPER-FILM POUCH

IMPORTANT INFO

STERILIZATION METHOD	GAMMA IRRADIATION
STORAGE TEMPERATURE	2/30°C – DO NOT FREEZE
SHELF LIFE	36 MONTHS

This bulletin (Rev.01-18/05/2020) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.