



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	179CS10
GENERAL DESCRIPTION	ANSE PER INOCULAZ.10 UL SOFT
DETAILED DESCRIPTION	1 outer box contain 4 inner boxes each composed of 50 paper-film pouches each containing 10 light blue 10 uL loops.

PRODUCT COMPONENTS

		SPECIFICATIONS
LOOP	LOOP 10UL SOFT LIGHT BLUE	TOTAL LENGTH: 201.00 mm \pm 1.00 mm EXTERNAL DIAMETER RING: 6.50 mm \pm 0.20 mm INTERNAL DIAMETER RING: 3.90 mm \pm 0.20 mm THICKNESS RING: 1.20 mm \pm 0.20 mm MATERIAL: POLYSTYRENE

PACKAGING SPECIFICATION

INNER BOX QTY	500
OUTER BOX QTY	2000
PACKAGING	PAPER-FILM POUCH

IMPORTANT INFO

STERILIZATION METHOD	GAMMA IRRADIATION
STORAGE TEMPERATURE	2/30°C – DO NOT FREEZE
SHELF LIFE	36 MONTHS

This bulletin (Rev.02-01/07/2020) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.