



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	349C
GENERAL DESCRIPTION	UTM-RT 1,5ML TUBE BULK
DETAILED DESCRIPTION	6 boxes containing 50 tube with 1,5 ml of UTM medium and red cap

PRODUCT COMPONENTS

MEDIUM		COMPOSITION
UTM MEDIUM		SUCROSE HBSS SOLUTION BOVINE SERUM ALBUMIN BUFFERED SOLUTION GELATIN AMINO ACIDS ANTIBIOTICS PHENOL RED
SAMPLING TUBE	CONICAL TUBE	SPECIFICATIONS
		TUBE LENGTH: 101.50 mm ± 0.20 mm TUBE DIAMETER: 16.30 mm ± 0.20 mm INSIDE DIAMETER: 14.10 mm ± 0.20 mm MATERIAL: POLYPROPYLENE
CAP	RED CAP	SPECIFICATIONS
		LOWER CAP DIAMETER: 21.50 mm ± 0.20 mm CAP LENGTH: 16.30 mm ± 0.20 mm HIGHER CAP DIAMETER: 20.50 mm ± 0.20 mm MATERIAL: HIGH DENSITY POLYETHYLENE

This bulletin (Rev.02-17/01/2020) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.



TECHNICAL DATA SHEET

PACKAGING SPECIFICATION

INNER BOX QTY	50
OUTER BOX QTY	300
PACKAGING	BULK

IMPORTANT INFO

STORAGE TEMPERATURE	2/25°C-DO NOT FREEZE
SHELF LIFE	18 MONTHS
MEDIUM pH	7.3 ± 0.2