



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	353C
GENERAL DESCRIPTION	UTM MINI 3ML IN BULK
DETAILED DESCRIPTION	6 boxes containing 50 tubes in bulk with 3 ml of UTM medium and red screw cap.

PRODUCT COMPONENTS

MEDIUM	COMPOSITION
UTM MEDIUM	SUCROSE HBSS SOLUTION BOVINE SERUM ALBUMIN BUFFERED SOLUTION GELATIN AMINO ACIDS ANTIBIOTICS PHENOL RED

SAMPLING TUBE	TUBE	SPECIFICATIONS
		TOTAL LENGTH: 81.30 mm ± 0.20 mm OUTER TUBE DIAMETER: 13.00 mm ± 0.20 mm INNER TUBE DIAMETER: 11.50 mm ± 0.20 mm MATERIAL: POLYPROPYLENE
CAP	RED CAP	SPECIFICATIONS
		CAP HEIGHT: 15.00 mm ± 0.20 mm LOWER CAP DIAMETER: 16.80 mm ± 0.20 mm HIGHER CAP DIAMETER: 16.00 mm ± 0.20 mm MATERIAL: HIGH DENSITY POLYETHYLENE

This bulletin (Rev.01-17/01/2020) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.

**TECHNICAL DATA SHEET****PACKAGING SPECIFICATION**

INNER BOX QTY	50
OUTER BOX QTY	300
PACKAGING	BULK

IMPORTANT INFO

STORAGE TEMPERATURE	5/25°C-DO NOT FREEZE
SHELF LIFE	18 MONTHS
MEDIUM pH	7.3 ± 0.2