



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	355C
GENERAL DESCRIPTION	UTM MINI 3ML + 1A/PL FLOCK REG
DETAILED DESCRIPTION	1 outer box contain 6 inner boxes each composed of 50 paper-film pouches each containing 1 tube with red screw cap filled with 3 ml of UTM medium and 1 flocked plastic regular swab individually packaged in a paper-film pouch.

PRODUCT COMPONENTS

MEDIUM	COMPOSITION
UTM MEDIUM	SUCROSE HBSS SOLUTION BOVINE SERUM ALBUMIN BUFFERED SOLUTION GELATIN AMINO ACIDS ANTIBIOTICS PHENOL RED

SAMPLING TUBE	TUBE	SPECIFICATIONS
		TOTAL LENGTH: 81.30 mm ± 0.20 mm OUTER TUBE DIAMETER: 13.00 mm ± 0.20 mm INNER TUBE DIAMETER: 11.50 mm ± 0.20 mm MATERIAL: POLYPROPYLENE
CAP	RED CAP	SPECIFICATIONS
		CAP HEIGHT: 15.00 mm ± 0.20 mm LOWER CAP DIAMETER: 16.80 mm ± 0.20 mm HIGHER CAP DIAMETER: 16.00 mm ± 0.20 mm MATERIAL: HIGH DENSITY POLYETHYLENE

This bulletin (Rev.03-13/05/2020) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.



TECHNICAL DATA SHEET

APPLICATOR	FLOCKED REGULAR SWAB	SPECIFICATIONS
		BREAK POINT: 77.50 mm \pm 1.00 mm STICK DIAMETER: 2.50 mm \pm 0.10 mm TOTAL LENGTH: 151.00 mm \pm 1.00 mm FLOCKED DIAMETER: 5.50 mm \pm 0.50 mm FLOCKED LENGTH: 16.00 mm \pm 2.00 mm STICK MATERIAL: POLYSTYRENE

PACKAGING SPECIFICATION

INNER BOX QTY	50
INNER BOX QTY	300
PACKAGING	PAPER-FILM POUCH

IMPORTANT INFO

STORAGE TEMPERATURE	2/25°C-DO NOT FREEZE
SHELF LIFE	18 MONTHS
MEDIUM pH	7.3 \pm 0.2

This bulletin (Rev.03-13/05/2020) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.