



## TECHNICAL DATA SHEET

### PRODUCT DESCRIPTION

<b>PRODUCT CODE</b>	359C
<b>GENERAL DESCRIPTION</b>	UTM MINI 1ML RT+ A/FLOC REG
<b>DETAILED DESCRIPTION</b>	6 boxes containing 50 paper-film pouch with one tube with red screw cap and 1 ml of UTM medium and 1 paper-film pouch containing one flocced plastic regular swab

### PRODUCT COMPONENTS

<b>MEDIUM</b>		<b>COMPOSITION</b>
UTM MEDIUM		SUCROSE HBSS SOLUTION BOVINE SERUM ALBUMIN BUFFERED SOLUTION GELATIN AMINO ACIDS ANTIBIOTICS PHENOL RED
<b>SAMPLING TUBE</b>	TUBE	<b>SPECIFICATIONS</b>
		TOTAL LENGTH: 81.30 mm ± 0.20 mm OUTER TUBE DIAMETER: 13.00 mm ± 0.20 mm INNER TUBE DIAMETER: 11.50 mm ± 0.20 mm MATERIAL: POLYPROPYLENE
<b>CAP</b>	RED CAP	<b>SPECIFICATIONS</b>
		CAP HEIGHT: 15.00 mm ± 0.20 mm LOWER CAP DIAMETER: 16.80 mm ± 0.20 mm HIGHER CAP DIAMETER: 16.00 mm ± 0.20 mm MATERIAL: HIGH DENSITY POLYETHYLENE

This bulletin (Rev.07-18/05/2020) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.



## TECHNICAL DATA SHEET

APPLICATOR	FLOCKED REGULAR SWAB	SPECIFICATIONS
		BREAK POINT: 77.50 mm $\pm$ 1.00 mm STICK DIAMETER: 2.50 mm $\pm$ 0.10 mm TOTAL LENGTH: 151.00 mm $\pm$ 1.00 mm FLOCKED DIAMETER: 5.50 mm $\pm$ 0.50 mm FLOCKED LENGTH: 16.00 mm $\pm$ 2.00 mm STICK MATERIAL: POLYSTYRENE

### PACKAGING SPECIFICATION

INNER BOX QTY	50
OUTER BOX QTY	300
PACKAGING	PAPER-FILM POUCH

### IMPORTANT INFO

STORAGE TEMPERATURE	2/25°C-DO NOT FREEZE
SHELF LIFE	18 MONTHS
MEDIUM pH	7.3 $\pm$ 0.2

This bulletin (Rev.07-18/05/2020) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.