



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	360C
GENERAL DESCRIPTION	UTM MINI 1ML RT+ A/FL.PERNASAL
DETAILED DESCRIPTION	6 boxes containing 50 paper film pouches each containing 1 tube with red cap filled with 1ml of UTM solution and 1 flocked swab packaged in paper-film pouch.

PRODUCT COMPONENTS

MEDIUM	COMPOSITION
UTM MEDIUM	SUGARS HBSS SOLUTION BOVINE SERUM ALBUMIN BUFFERED SOLUTION JELLY AMINO ACIDS ANTIBIOTICS PH INDICATOR

SAMPLING TUBE	TUBE	SPECIFICATIONS
		TOTAL LENGTH: 81.30 mm ± 0.20 mm TUBE DIAMETER: 13.00 mm ± 0.20 mm INTERNAL DIAMETER: 11.50 mm ± 0.20 mm MATERIAL: POLYPROPYLENE
CAP	RED CAP	SPECIFICATIONS
		CAP HEIGHT: 15.00 mm ± 0.20 mm LOWER CAP DIAMETER: 16.80 mm ± 0.20 mm HIGHER CAP DIAMETER: 16.00 mm ± 0.20 mm MATERIAL: HIGH DENSITY POLYETHYLENE

This bulletin (Rev.05-26/08/2019) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.



TECHNICAL DATA SHEET

APPLICATOR	FLOCKED SWAB	SPECIFICATIONS
		BREAKING POINT: 98.50 mm \pm 1.00 mm DIAMETER HANDLE PART: 2.50 mm \pm 0.10 mm TOTAL LENGTH: 152.00 mm \pm 2.00 mm FLOCKED DIAMETER: 4.00 mm \pm 0.50 mm FLOCKED LENGTH: 16.00 mm \pm 2.00 mm STICK MATERIAL: POLYSTYRENE

PACKAGING SPECIFICATION

INNER BOX QTY	50
OUTER BOX QTY	300
PACKAGING	PAPER-FILM POUCH

IMPORTANT INFO

STORAGE TEMPERATURE	2/25°C-DO NOT FREEZE
SHELF LIFE	18 MONTHS
MEDIUM pH	7.3 \pm 0.2