



## TECHNICAL DATA SHEET

### PRODUCT DESCRIPTION

<b>PRODUCT CODE</b>	365C
<b>GENERAL DESCRIPTION</b>	UTM MINI 2ML PERNASAL FLOCKED
<b>DETAILED DESCRIPTION</b>	6 boxes containing 50 bags paper-film each containing 1 tube with 2 ml UTM medium and red cap, 1 flocked pernasal swab in paper-film pouch

### PRODUCT COMPONENTS

<b>MEDIUM</b>	<b>COMPOSITION</b>
UTM MEDIUM	SUCROSE HBSS SOLUTION BOVINE SERUM ALBUMIN BUFFERED SOLUTION GELATIN AMINO ACIDS ANTIBIOTICS PHENOL RED

<b>SAMPLING TUBE</b>	<b>TUBE</b>	<b>SPECIFICATIONS</b>
		TOTAL LENGTH: 81.30 mm ± 0.20 mm OUTER TUBE DIAMETER: 13.00 mm ± 0.20 mm INNER TUBE DIAMETER: 11.50 mm ± 0.20 mm MATERIAL: POLYPROPYLENE
<b>CAP</b>	<b>RED CAP</b>	<b>SPECIFICATIONS</b>
		CAP HEIGHT: 15.00 mm ± 0.20 mm LOWER CAP DIAMETER: 16.80 mm ± 0.20 mm HIGHER CAP DIAMETER: 16.00 mm ± 0.20 mm MATERIAL: HIGH DENSITY POLYETHYLENE

This bulletin (Rev.02-17/01/2020) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.



## TECHNICAL DATA SHEET

APPLICATOR	FLOCKED SWAB	SPECIFICATIONS
		BREAK POINT: 98.50 mm $\pm$ 1.00 mm STICK DIAMETER: 2.50 mm $\pm$ 0.10 mm TOTAL LENGTH: 152.00 mm $\pm$ 2.00 mm FLOCKED DIAMETER: 4.00 mm $\pm$ 0.50 mm FLOCKED LENGTH: 16.00 mm $\pm$ 2.00 mm STICK MATERIAL: POLYSTYRENE

### PACKAGING SPECIFICATION

INNER BOX QTY	50
OUTER BOX QTY	300
PACKAGING	PAPER-FILM POUCH

### IMPORTANT INFO

STORAGE TEMPERATURE	2/25°C-DO NOT FREEZE
SHELF LIFE	18 MONTHS
MEDIUM pH	7.3 $\pm$ 0.2