



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	366C
GENERAL DESCRIPTION	UTM MINI 2ML REGULAR FL.APP.
DETAILED DESCRIPTION	6 boxes containing 50 paper-film pouches containing 1 flocked plastic regular swab in a paper-film pouch and a tube with red screw cap and 2 ml of UTM medium

PRODUCT COMPONENTS

MEDIUM	COMPOSITION
UTM MEDIUM	SUCROSE HBSS SOLUTION BOVINE SERUM ALBUMIN BUFFERED SOLUTION GELATIN AMINO ACIDS ANTIBIOTICS PHENOL RED

SAMPLING TUBE	TUBE	SPECIFICATIONS
		TOTAL LENGTH: 81.30 mm ± 0.20 mm OUTER TUBE DIAMETER: 13.00 mm ± 0.20 mm INNER TUBE DIAMETER: 11.50 mm ± 0.20 mm MATERIAL: POLYPROPYLENE
CAP	RED CAP	SPECIFICATIONS
		CAP HEIGHT: 15.00 mm ± 0.20 mm LOWER CAP DIAMETER: 16.80 mm ± 0.20 mm HIGHER CAP DIAMETER: 16.00 mm ± 0.20 mm MATERIAL: HIGH DENSITY POLYETHYLENE

This bulletin (Rev.02-17/01/2020) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.



TECHNICAL DATA SHEET

APPLICATOR	FLOCKED REGULAR SWAB	SPECIFICATIONS
		STICK DIAMETER: 2.50 mm \pm 0.10 mm TOTAL LENGTH: 151.00 mm \pm 2.00 mm FLOCKED DIAMETER: 5.50 mm \pm 0.50 mm FLOCKED LENGTH: 16.00 mm \pm 2.00 mm BREAK POINT: 78.50 mm \pm 1.00 mm STICK MATERIAL: POLYSTYRENE

PACKAGING SPECIFICATION

INNER BOX QTY	50
OUTER BOX QTY	300
PACKAGING	PAPER-FILM POUCH

IMPORTANT INFO

STORAGE TEMPERATURE	2/25°C-DO NOT FREEZE
SHELF LIFE	18 MONTHS
MEDIUM pH	7.3 \pm 0.2

This bulletin (Rev.02-17/01/2020) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.