



## TECHNICAL DATA SHEET

### PRODUCT DESCRIPTION

<b>PRODUCT CODE</b>	503CS01
<b>GENERAL DESCRIPTION</b>	FLEXIBLE FLOQSWAB SINGLE WRAPPED BP100
<b>DETAILED DESCRIPTION</b>	1 outer box contain 10 inner boxes each composed of 100 flexible flocked swabs individually packaged in paper-film pouch.

### PRODUCT COMPONENTS

		<b>SPECIFICATIONS</b>
<b>APPLICATOR</b>	FLEXIBLE FLOCKED SWAB	TOTAL LENGTH: 152.00 mm ± 2.00 mm STICK DIAMETER: 2.50 mm ± 0.10 mm FLOCKED LENGTH: 16.00 mm ± 2.00 mm FLOCKED DIAMETER: 4.00 mm ± 0.50 mm BREAK POINT: 98.50 mm ± 1.00 mm MATERIAL: POLYSTYRENE

### PACKAGING SPECIFICATION

<b>INNER BOX QTY</b>	100
<b>OUTER BOX QTY</b>	1000
<b>PACKAGING</b>	PAPER-FILM POUCH

### IMPORTANT INFO

<b>STERILIZATION METHOD</b>	ETO TREATMENT
<b>STORAGE TEMPERATURE</b>	2/30°C – DO NOT FREEZE
<b>SHELF LIFE</b>	36 MONTHS

This bulletin (Rev.13-18/08/2020) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.