



TECHNICAL DATA SHEET

PRODUCT DESCRIPTION

PRODUCT CODE	516CS01
GENERAL DESCRIPTION	MINI FLEXIBLE FLOQSWAB SIGLE WRAPPED BP100
DETAILED DESCRIPTION	1 outer box contain 10 inner boxes each composed of 100 flocked swabs individually wrapped in paper-film pouch.

PRODUCT COMPONENTS

		SPECIFICATIONS
APPLICATOR	FLEXIBLE FLOCKED SWAB	TOTAL LENGTH: 152.00 mm ± 2.00 mm STICK DIAMETER: 2.50 mm ± 0.10 mm BREAKING POINT: 98.50 mm ± 1.00 mm FLOCKED LENGTH: 8.00 mm ± 2.00 mm FLOCKED DIAMETER: 3.50 mm ± 0.50 mm STICK MATERIAL: POLYSTYRENE

PACKAGING SPECIFICATION

INNER BOX QTY	100
OUTER BOX QTY	1000
PACKAGING	PAPER-FILM POUCH

IMPORTANT INFO

STERILIZATION METHOD	ETO TREATMENT
STORAGE TEMPERATURE	2/30°C – DO NOT FREEZE
SHELF LIFE	36 MONTHS

This bulletin (Rev.07-16/11/2020) is for your guidance and is based on information believed to be reliable and should be periodically reviewed.